



KNITTED WRIST SPLINT

OS 3200

It is manufactured, in compliance with the wrist anatomy, from breathable elastic knitted fabric. There is a plastic splint supporting on the volar surface.

Indication

- Carpal tunnel syndrome (CTS)
- In Rheumatic arthritis
- For the purposes of stabilizing the wrist

Size

S	M	L	XL	XXL
14-15	15-16	16-17	17-18	18-19

< LEFT RIGHT >



OS 4300

KNITTED LUMBOSTAD CORSET

It is manufactured from intensively air permeable fabric made of cotton and knitted fabric. It is supported by plastic and steel stays in the lumbar and sacral regions.



Indication

- Lumbago
- Spinal deformities
- Vertebral and ilio-sacral joint problems

Size

S	M	L	XL	XXL
80-90	90-100	100-110	110-120	120-130



KNITTED LUMBOSACRAL CORSET OS 4300 A

It is manufactured, in compliance with the lumbar region anatomy, from breathable elastic knitted fabric. It is supported by plastic and steel stays in the lumbar and sacral regions.

Indication

- Lumbago
- Spinal deformities
- Vertebral and ilio-sacral joint problems

Size

S	M	L	XL	XXL
80-90	90-100	100-110	110-120	120-130





KNITTED ANKLE GUARD WITHOUT SUPPORT

OS 1514

It is manufactured, in compliance with the ankle anatomy, from breathable elastic knitted fabric. It is manufactured in order to stabilize the ankle.

Indication

- For thermo therapy purposes
- In cases where stabilization is required
- In the dissipation of post-distortion edema



Size

S	M	L
18-20	22-24	26-28

OS 1500

KNITTED ANKLE GUARD WITH MALLEOLUS SUPPORT

It is manufactured, in compliance with the ankle anatomy, from breathable elastic knitted fabric. There is a silicone pad supporting medial and lateral malleoli.



Indication

- In the dissipation of post-distortion edema
- For thermo therapy purposes

Size

S	M	L	XL	XXL
18	20	22	24	26



KNITTED ANKLE GUARD WITH ACHILLES TENDON SUPPORT

OS 1518

It is manufactured, in compliance with the ankle anatomy, from breathable elastic knitted fabric. There are a silicone pads supporting Achilles tendon and heel region.

Indication

- Achilles tendon injuries
- Post-Achilles tendon operation
- Achilles tendon inflammation



Size

S	M	L	XL	XXL
18	20	22	24	26





KNITTED KNEE BRACE WITH OUT SUPPORT

OS 1416

Indication

- For Thermo therapy purpose
- In cases where stabilization is required.
- In the dissipation of post-distortion edema
- Arthritis



Size

S	M	L
28-31	31-33	33-36

OS 1400

KNITTED KNEE BRACE WITH LIGAMENT SUPPORT

It is a type having silicone support for patella and ligament and manufactured from knitted material and is an air permeable anti allergen knee brace.



Indication

- Arthritis
- Effusion
- When the knee is painful
- For Thermo therapy purposes

Size

S	M	L	XL	XXL
28	30	32	34	36



KNITTED ELBOW BRACE WITH EPICONDYLITIS SUPPORT

OS 1216

It is manufactured, in compliance with the elbow region anatomy, from breathable elastic knitted fabric, It is supported with silicone pad from medial and lateral regions.

Indication

- Tenosynovitis
- Arthralgia
- In case of pain in the elbow
- In condylar injuries



Size

S	M	L	XL	XXL
19-21	21-23	23-25	25-27	27-29

