



## LATEX BANDI

05 9130

It is made of latex material with different hardness for exercise and rehabilitation purposes. Available for sale in 45.5 m in rolls or 1.20 cm as packs.

### Indication

- For exercise and
- Rehabilitation purposes

### Color



05 9172


## SILICONE EXERCISE BALL

It is manufactured, in compliance with the hand anatomy, in three different resistances from silicone

### Indication

- For hand exercise purposes
- For Rehabilitation purposes

### Hardness



soft

medium

hard



## HAND EXERCISE BALL (SPONGE)

05 9009

Made from sponge material in compliance with the hand anatomy.

### Indication

- For hand exercise purposes
- For Rehabilitation purposes

### Size

STANDARD



## EPIN SUPPORTED CORK INSOLE

**OS 5509**

It is manufactured, in compliance with the foot ergonomics, offering longitudinal arch support with metatarsal support using epine cushioned cork and artificial leather.

### Indication

- Fallen arch
- Pes Planus (flat-footedness)
- Pes plano valgus



### Size

34 35 36 37 38 39 40 41 42 43 44 45 46

**OS 5503**

## SILICONE HALF INSOLE 3/4

It is Made of 100% silicone. It is Hypo-allergenic. Blue areas are softer than other areas. It extends up to the metatarsal region.



### Indication

- In case of the following conditions •
- Plantar fasciitis
  - Metatarsalgia
  - Pes planus
- As a shock absorber •  
Keeping Achilles tendon intact

### Size

S	M	L	XL
34-35-36	37-38-39	40-41-42	43-44-45



## SILICONE INSOLE

**OS 5504**

It is Made of 100% silicone. It is Hypo-allergenic. Thanks to the more elastic nature of the blue areas than the other areas it keeps the muscles in the foot sole, active.

### Indication

- In case of the following conditions
- Plantar fasciitis
- Metatarsalgia
- Pes planus
- As a shock absorber
- Keeping achilles tendon intact



### Size

XS	S	M	L	XL	XXL
31-33	34-36	37-39	40-42	43-45	46-48





## SILICONE INSOLE WITH INTENSIVE ARCHSUPPORT

**OS 5514**

It is Made of 100% silicone. It is Hypo-allergenic. Thanks to the more elastic nature of the blue areas than the other areas it keeps the muscles in the foot sole, active.

### Indication

- Plantar faciitis
- Metatarsalgia
- Pes planus
- As a shock absorber
- Keeping Achilles tendon intact



### Size

XS	S	M	L	XL	XXL
35-36	37-38	39-40	41-42	43-44	45-46

**OS 5505**

## STEELINSOLE

It is used in the treatment of Pes Planus, metatarsalgia and deformities. It is provided with a height-adjustable stainless steel bridge and distinguished supports. It maintains its elasticity and dynamics for a long time. it was made from the natural cowhide leather.



### Size

34	35	36	37	38	39	40	41	42	43	44	45	46
----	----	----	----	----	----	----	----	----	----	----	----	----

## ORTHOFLEX INSOLE

**OS 5506**

It is manufactured, in compliance with the foot ergonomics, from epine cushioned plastazote artificial leather with longitudinal arch and metatarsal support.

### Indication


- Fallen arch
- Pes planus
- Epine calcaneal
- Pes plano valgus



### Size

34	35	36	37	38	39	40	41	42	43	44	45	46
----	----	----	----	----	----	----	----	----	----	----	----	----





## SILICONE HEEL CUP

**OS 5502**

It is manufactured, in compliance with the foot anatomy, from 100% silicone. It is Hypo-allergenic. blue areas are softer than the other areas.

### Indication

- Heel spur
- General heel pain
- Plantar fasciitis
- irritations in the Achilles tendon



### Size

S	M	L
35-38	39-42	43-44

## OS 9004

## UCBL INSOLE

UCBL, which is used in children patients, distributes the whole weight of the body equally to the feet, with its rigid structural frame and unique design. It grips the heel bone and supports in a fully orthopedic manner with its proper design.



### Indication

- In case of the following conditions •
- Plantar fasciitis
  - Metatarsalgia
  - Pes planus
  - As a shock absorber
  - Keeping Achilles tendon intact

### Size

XS	S	M	L	XL	XXL
24-25	26-27	28-29	30-31	32-33	34-35



## SILICONE MEDIAL/LATERAL WEDGE (ADHESIVE)

**OS 9059**

It is manufactured, in compliance with the foot anatomy, from silicone coated with fabric, It has a special adhesive at the bottom surface to stick to the shoe. It can be used in both directions.

### Indication

- It is used in order to prevent from in-toeing and out-toeing.



### Size

S	M	L	XL
34-35-36	37-38-39	40-41-42	43-44-45





## SILICONE REEL TOE SEPARATOR

**OS 5501**

It is manufactured from 100% silicone.  
It is Hypo-allergenic. Available in blue and white color.

### Indication

- Positioning for Hallux valgus



### Size

<b>S</b>	<b>M</b>	<b>L</b>
35-38	39-42	43-44

**OS 9073 A**

## PAD TOE SEPARATOR FOR BUNIONS

It is manufactured from 100% silicone. It is Hypo-allergenic. It has a reel to separate the big toe and pad to support bunion.



### Indication

- Positioning for Hallux valgus • protection of bunion



### Size

Big Toe      Little Toe



## SILICONE HEEL SOCK

**OS 5519**

It is manufactured, in compliance with the heel anatomy, from 100% silicone. It is Hypo-allergenic. It provides comfort by wearing it on heel like a sock.

### Indication

- Heel spur
- Comfort



### Size

STANDARD



## ARCH SUPPORT WITH REMOVABLE ADHESIVE

**OS 5522**

It is manufactured from fabric-coated silicone to support the longitudinal arch region. There is a special adhesive underneath the silicone for sticking the shoe.

### Indication

- Pes planus
- Supporting arch longitudinally



### Size

<b>S</b>	<b>M</b>	<b>L</b>
34-37	38-41	41-45

**OS 5523**

## SELF-ADHESIVE EPINPAD

It is manufactured from fabric-coated silicone to support the heel region. There is a special adhesive underneath the silicone for sticking the shoe.



### Indication

- Heel spur
- General heel pain
- Plantar fasciitis
- Irritations in the Achilles tendon

### Size

<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
34-35-36	37-38-39	40-41-42	43-44-45



## ADHESIVE METATARSAL PAD

**OS 5524**

It is manufactured from fabric-coated silicone to support the metatarsal region. There is a special adhesive underneath the silicone for sticking the shoe.

### Indication

- Metatarsal head sensitivities
- Metatarsal head callosities
- Metatarsalgia



### Size

<b>S</b>	<b>M</b>	<b>L</b>
34-37	38-41	41-45





## METATARSAL PAD WITH RING

**OS 9074**

In case of supporting and protecting the metatarsal area.  
There is a ring to fix the foot.

### Indication

- Metatarsal head sensitivities
- Metatarsal head callosities
- Metatarsalgia



### Size

STANDARD

**OS 9075**

## METATARSAL PAD WITH SEPARATOR

In case of supporting and protecting the metatarsal area.  
There is a wedge to fix the foot.



### Indication

- Metatarsal head sensitivities
- Metatarsal head callosities
- Metatarsalgia



### Size

STANDARD



## EPINE SPONGE HEEL CUP

**OS 9028**

It is manufactured from fabric-coated foam to support the heel region. There is a special adhesive underneath the silicone for sticking the shoe. there is a removable part in the middle section.

### Indication

- Heel spur



### Size

S

M

34-39

40-45