



## PROPHYLACTIC COMPRESSION KNEE-HIGH STOCKING (PROFILEACTIVE)

**OS 5001**

Orthosoft compression stockings are designed and manufactured in accordance with TS ENV 12718 standard. Graduated pressure feature make contribution to your health by treating your legs. Special and high quality raw materials are used in the production. It has a soft touch around it, it is easy to wear and take off.


AD  
Knee-high

Pressure Class



140 Denier  
Profile Active  
8-12 mmHg  
pressure

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



**OS 5002**

## PROPHYLACTIC COMPRESSION THIGH-HIGH STOCKING (PROFILE ACTIVE)

Orthosoft compression stockings are designed and manufactured in accordance with TS ENV 12718 standard. Graduated pressure feature make contribution to your health by treating your legs. Special and high quality raw materials are used in the production. It has a soft touch around it, it is easy to wear and take off.


AG  
Thigh-high

Pressure Class



140 Denier  
Profile Active  
8-12 mmHg  
pressure

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



## PROPHYLACTIC COMPRESSION PANTY HOSE (PROFILE ACTIVE)

**OS 5003**

Orthosoft compression stockings are designed and manufactured in accordance with TS ENV 12718 standard. Graduated pressure feature make contribution to your health by treating your legs. Special and high quality raw materials are used in the production. It has a soft touch around it, it is easy to wear and take off.


AT  
PantyHose

Pressure Class



140 Denier  
Profile Active  
8-12 mmHg  
pressure

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79





## ANTIEMBOLI KNEE-HIGH STOCKINGS

**OS 5006**

Orthosoft Anti-Embolism stockings are designed and produced according to TS ENV 12719 standard. It offers a graduated pressure starting from the ankle to the thigh.

Special and high quality raw materials are used in production. Each type of stocking has a control hole.


AD  
Knee-high

Pressure Class



Pressure at the anti-embolism ankle 18-21 mmHg.

Lower Calf	Size				
	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32
cD	26-32	30-36	34-40	37-43	40-47
cG	43-57	49-63	53-67	58-73	64-79



**OS 5005**

## ANTIEMBOLI THIGH-HIGH STOCKINGS

Orthosoft Anti-Embolism stockings are designed according to TSENV 12719 standard. It offers a graduated pressure starting from the ankle to the thigh.

Special and high quality raw materials are used in production. Each type of stocking has a control hole.

At the upper section, special drop silicone rubbers are used which prevent sliding over the skin and introducing a turnstile


AG  
Thigh-high

Pressure Class



Pressure at the anti-embolism ankle 18-21

Lower Calf	Size				
	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32
cD	26-32	30-36	34-40	37-43	40-47
cG	43-57	49-63	53-67	58-73	64-79



## ANTIEMBOLI COMPRESSION STOCKINGS WITH WAIST ATTACHMENT

**OS 5007**

Orthosoft Anti-Embolism stockings are designed according to TSENV 12719 standard. It offers a graduated pressure starting from the ankle to the thigh.

Special and high quality raw materials are used in production. Each type of stocking has a control hole. Those with waist attachment at the waist have adjustable belts.


AT  
Belted waist

Pressure Class



Pressure at the anti-embolism ankle 18-21 mmHg.

Lower Calf	Size				
	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32
cD	26-32	30-36	34-40	37-43	40-47
cG	43-57	49-63	53-67	58-73	64-79





## CCL 1 KNEE-HIGH COMPRESSION STOCKINGS (LOW PRESSURE)

**OS 5010**

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it, it is easy to wear and take off.

✓ OPEN ✓ CLOSED


AD  
Knee-high

Pressure Class



CCL1 Low pressure at the ankle 18-21 mmHg.

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



**OS 5011**

## CCL 1 KNEE-HIGH COMPRESSION STOCKINGS (LOW PRESSURE)

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it, it is easy to wear and take off.

✓ OPEN ✓ CLOSED


AG  
Thigh-high

Pressure Class



CCL1 Low pressure at the ankle 18-21 mmHg.

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



## CCL 1 PANTYHOSE COMPRESSION STOCKINGS (LOW PRESSURE)

**OS 5012**

Orthosoft Anti-Embolism stockings are designed according to TSENV 12719 standard. it offers a graduated pressure starting from the ankle to the thigh. Special and high quality raw materials are used in production. Each type of stocking has a control hole.

✓ OPEN ✓ CLOSED


AT  
Pantyhose

Pressure Class



CCL1 Low pressure at the ankle 18-21 mmHg.

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79





## CCL2 KNEE-HIGH COMPRESSION STOCKINGS (MEDIUM PRESSURE)

**OS 5013**

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it, it is easy to wear and take off.

✓ OPEN ✓ CLOSED


AD  
Knee-high

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



**OS 5014**

## CCL2 THIGH-HIGH COMPRESSION STOCKINGS (MEDIUM PRESSURE)

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it. it is easy to wear and take off.

✓ OPEN ✓ CLOSED


AG  
Thigh-high

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



## CCL2 PANTY HOSE COMPRESSION STOCKINGS (MEDIUM PRESSURE)

**OS 5015**

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it, it is easy to wear and take off.

✓ OPEN ✓ CLOSED


AT  
Panty hose

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79





## PANTY HOSE COMPRESSION STOCKING (FOR PREGNANT)

**OS 5009**

Orthosoft compression stockings are designed and produced according to TS ENV 12718 standard. Thanks to graduated pressure, it makes contribution to your health by treating your legs. Special and high quality raw materials are used in production. It has a soft touch around it, it is easy to wear and take off.

✓ OPEN ✓ CLOSED


ATU pregnant panty

Pressure Class



CCL1 Low pressure at the ankle 18-21 mmHg.

Lower Calf	Size						
	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cD	27-32	28-35	32-38	34-41	36-44	38-47	40-50
cG	43-58	46-61	49-64	52-68	56-71	60-74	64-79



**OS 5016**

## SILICONE COMPRESSION ARMSLEEVE WRIST-SHOULDER


Lymph edema Compression arm sleeves are designed and produced according to RAL GZ 387/2 standard. The pressure exerted recovers the damaged lymphatic vessels and removal of the lymph nodes in chest surgery prevents the clogging in the circulation afterwards. Special and high quality materials are used in the production. It has a soft touch. it is easy to wear and take off.



Silicone Pressure Class

CCL2 medium pressure at the ankle 23-32 mmHg

Arm circumference	Size						
	1	2	3	4	5	6	7
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cE	16-22	21-27	24-30	27-33	30-36	33-39	36-42
cG	20-27	23-30	26-33	30-37	34-41	38-45	42-49



## SILICONE COMPRESSION ARMSLEEVE FINGER-SHOULDER

**OS 5017**

Compression sleeves are designed and manufactured in accordance with RAL GZ 387/2. The pressure exerted recovers the damaged lymphatic vessels and removal of the lymph nodes in chest surgery prevents the clogging in the circulation afterwards. Special and high quality materials are used in the production. It has a soft touch. it is easy to wear and take off.


Silicone

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Arm circumference	Size						
	1	2	3	4	5	6	7
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cE	16-22	21-27	24-30	27-33	30-36	33-39	36-42
cG	20-27	23-30	26-33	30-37	34-41	38-45	42-49





## COMPRESSION ARMSLEEVE WITH SHOULDER STRAP / WRIST-SHOULDER

**OS 5018**

Lymph edema Compression arm sleeves are designed and produced according to RAL GZ 387/2 standard. The pressure exerted recovers the damaged lymphatic vessels and removal of the lymph nodes in chest surgery prevents the clogging in the circulation afterwards. Special and high quality materials are used in the production. It has a soft touch. it is easy to wear and take off.


Silicone

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Arm circumference	Size						
	1	2	3	4	5	6	7
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cE	16-22	21-27	24-30	27-33	30-36	33-39	36-42
cG	20-27	23-30	26-33	30-37	34-41	38-45	42-49



**OS 5019**

## COMPRESSION ARMSLEEVE W/ STRAP FINGER-SHOULDER

Lymph edema Compression arm sleeves are designed and produced according to RAL GZ 387/2 standard. The pressure exerted recovers the damaged lymphatic vessels and removal of the lymph nodes in chest surgery prevents the clogging in the circulation afterwards. Special and high quality materials are used in the production. It has a soft touch. it is easy to wear and take off.


Silicone

Pressure Class



CCL2 medium pressure at the ankle 23-32 mmHg

Arm circumference	Size						
	1	2	3	4	5	6	7
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cE	16-22	21-27	24-30	27-33	30-36	33-39	36-42
cG	20-27	23-30	26-33	30-37	34-41	38-45	42-49



## COMPRESSION STOCKING WEARING APPARATUS

**OS 9228**

It is difficult to wear since the varicose stockings are pressurized. Thanks to this product, socks can be worn easily without difficulty in a very short time.

It is manufactured from a material that will not harm metal and painted socks. It provides comfort in wearing compression stockings.