



## WRIST SPLINT

GR 1201 Ç

It is manufactured, in compliance with the child anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There is an aluminum splint supporting on the volar surface. no distinction such as right or left.

### Indication

- Carpal tunnel syndrome (CTS)
- Rheumatic arthritis
- Painful condition of the wrist



### Size

STANDARD

GR 1202 Ç

## WRIST SPLINT WITH THUMB SUPPORT

It is manufactured, in compliance with the child anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There is an aluminum splint supporting on the volar surface. No distinction such as right or left.



### Indication

- De Quarvein
- Carpal tunnel syndrome (CTS)
- Rheumatic arthritis
- Painful condition of the wrist

### Size

STANDARD



## DE QUVARIN SPLINT

GR 1203 Ç

It is manufactured, in compliance with the child anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There is an aluminum splint supporting on the thenar region.

### Indication

- De Quarvein
- 1 .CMC and 1 .MCP traumas
- In rheumatic conditions



### Size

STANDARD



## PATELLASTABILIZING KNEE BRACE

GR 1403 Ç

It is manufactured, in compliance with the child anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There is a circular pad around patella for supporting purposes.

### Indication

- In case of Arthritis
- Effusion
- In case of painful knee
- For thermotherapy purposes



### Size

STANDARD

OS 1407 Ç

## ANKLE GUARD WITH LIGAMENT SUPPORT

It is manufactured, in compliance with the child knee anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There are flexible aluminum springs in order to support the ligaments in the mediolateral region.



### Indication

- In cases where compression pad application is required
- In cases where stabilization is required
- In distortion
- For Thermotherapy purposes

### Size

STANDARD



## LUMBOSTAD CORSET

OS 4304 Ç

It is manufactured, in compliance with the child anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. It has a mediolateral supporting pad. Its height is 27 cm. It is supported with plastic and steel stays in the lumbar region. extra belts are used in the waist region for better fit.

### Indication

- Lumbago
- Spinal deformities
- Vertebral problems

### Size

STANDARD





## ANKLE SUPPORT WITH 8 STRAPS **GR 1502 Ç**

It is manufactured, in compliance with the child ankle anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. Ankle posture is adjusted with extra elastic strap.

### Indication

- In the situations requiring stabilization
- In the dissipation of post-distortion edema
- For thermotherapy purposes



### Size

STANDARD

## **GR 1503 Ç**

## ANKLE GUARD WITH LIGAMENT SUPPORT

It is manufactured, in compliance with the child ankle anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. There are flexible aluminum springs in order to support the ligaments in the mediolateral region



### Indication

- Where compression pad application is required
- Where stabilization is required
- In case of distortion
- For Thermotherapy purposes

### Size

STANDARD



## ANKLE GUARD WITH PLASTIC SUPPORT **GR 1504 Ç**

It is manufactured, in compliance with the child ankle anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton. plastic splint exposing the malleoli is added in order to support ligaments.

### Indication

- In mediolateral ligament injuries
- In cases where Stabilization is required
- In dissipation of post-distortion edema
- For thermotherapy purposes



### Size

STANDARD



## ANKLE GUARD WITH MALLEOLUS SUPPORT

**GR 1505 Ç**

It is manufactured, in compliance with the child ankle anatomy, from 2.5 mm adhesive raised neoprene velcro fabric the inner surface of which is cotton fabric. The pads located around malleoli help dissipate the edema while moving.

### Indication

- Thendo Myopathies
- In distortion
- For thermotherapy purposes
- In cases where Stabilization is required



### Size

STANDARD

**OS 9039 Ç**

## CHILD TRUSS

It is employed in hernia cases before and after surgical operation. It prevents the growth of the hernia before surgery and the bursting of the stitches after surgery. It does not harm the body owing to its specially designed pillow.



### Size

XS-S	M-L	XL-XXL
80-90	95-110	115-125



## MINERVA ORTHOSIS (CHILD SIZE)

**OS 9027 Ç**

It is manufactured, in compliance with the child neck and body anatomy, from polyethylene material the inner surface of which is coated with fabric. It provides a rigid position to stabilize the cervical region. The height can be adjusted. It will not allow rotation.

### Indication

- In case of,
- Cervical region fractures,
- Cervical region trauma,
- Cervical region muscle and ligament injuries. and
- Cervical joint diseases.



### Size

STANDARD





## MESH ARM SLING (CHILD SIZE)

OS 2206 Ç

It is manufactured from special perforated fabric. There is a strap that allows the arm to stay at the desired position. It is easy and comfortable to use with its special buckle.

### Indication

- When the arm needs to be suspended
- In soft tissue injuries
- When the arm needs to remain stationary



### Size

XXS

XS

20-25

25-30

OS 2201 Ç

## SPONGE ARM SLING (CHILD SIZE)

It is manufactured from the sponge laminated with fabric the inner surface of which is cotton towel. Adhesive velcro fabric is used at the outer surface.



### Indication

For resting purposes in forearm problems •

### Size

XXS

XS

20-25

25-30



## VELPAUBANDAGE (CHILD SIZE)

OS 2208 Ç

It is manufactured from sponge laminated with fabric, the inner surface of which is cotton towel. Extra bandage is used to stabilize the shoulder joint.

### Indication

- In Shoulder- luxation, fracture and subluxation
- For resting purposes in case of forearm problems



### Size

STANDARD



## CLAVICLE BANDAGE (CHILD SIZE)

OS 2205 Ç

It is manufactured from the sponge laminated with fabric the inner surface of which is cotton towel.

### Indication

- For positioning epiphyses in clavicle fracture
- For stabilization purposes



### Size

XXS

XS

30-35

35-40

OS 4308 Ç

## POSTUREX CORSET (CHILD SIZE)

It is manufactured as 26 cm from intensively air permeable fabric made of cotton and fishing line. It is supported by plastic stays for dorsal adaptation.



### Indication

- In vortal
- In Postural disorders
- muscle weakness

### Size

XXS

XS

67-73

74-81



## CHILD WAIST BELT

OS 9107 Ç

It is manufactured as 6 cm in width from intensively air permeable fabric made of cotton and fishing line. There is a hard chock to put pressure on the belly.

### Indication

- Umbilical hernia



### Size

STANDARD





## FREJKA PILLOW (HIP ABDUCTION CORSET)

**OS 9044 Ç**

It is made of cotton fabric, and the sections in contact with the legs are supported with sponge. It allows the baby to do the compulsory hip and leg movements because a flexible material is used. It provides convenience while wearing with 4 pieces of adjustable step plastic buckles.

### Indication

- In case of hip dislocation

### Size

<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XL</b>
18	22	24	26	28	30

**OS 9041 Ç**

## PAVLIK BANDAGE

It is manufactured from velcro fabric with toweling fabric inside. Presence of adhesive velcro fabric at the outer surface ensures full application to the body. The buckles in addition to their stylish look, also make the such adjustment easy.

### Indication

- In case of hip dislocation •

### Size

<b>XS-S</b>	<b>M-L</b>	<b>XL-XXL</b>
35-45	45-55	55-65



## HIP ABDUCTION SHORTS (PLASTAZOTE)

**OS 9024 Ç**

It is manufactured, in compliance with the child hip anatomy, from plastazote material and immobilizes the hip at an adduction angle of 180° and a flexion angle of 30°.

### Indication

- In case of hip dislocation

### Size

<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
----------	----------	----------	----------	----------



## DENNIS BRAUN SPLINT

05 9022

It can allow the movements where the legs are independent of each other, such as crawl. It prevents the boot from being disengaged, and avoids also the regional rednesses that will occur on the heel.

### Indication

- In case of hip dislocation



### Size

S	M
29-41	42-68

05 9022 A

## DENNIS BROUNSHOE

It is a shoe designed to make use of Dennis Braun splint. It is manufactured from artificial leather. Its instep section is open. the leg section is high to grip the ankle comfortably.



### Indication

- In case of hip dislocation



### Size

12	13	14	15	16	17	18	19
27	26	25	24	23	22	21	20



## HIP ABDUCTION ORTHOSIS (PLASTIC)

05 9143

It is manufactured as a single piece from polyurethane material there is a cotton fabric inside. Its use is easy due to velcro. It immobilizes the hip at an adduction angle of 180° and a flexion angle of 30°.

### Indication

- Hip dislocation

### Size / Moon

S	M	L
0-2	3-5	6-10

05 9209

## HIP DYSPLASIA ORTHOSIS



### Indication

.....•

.....•

### Size / Moon

S	M
0-8	8-16







## CIRCUMCISION PANTIES

OS 9166 Ç

It is manufactured from cotton fabric and plastic material. It is used after circumcission to prevent the child's penis from sticking to the skin and from being damaged.

### Indication

- After circumcission

### Size

S

M

L

XL

XXL

## OS 9004

## UCBL INSOLES

UCL, which is used in children patients, distributes the whole weight of the body equally to the feet, with its rigid structural frame and unique design. It grips the heel bone and supports in a fully orthopedic manner with its proper design.



### Indication

- Plantar fasciitis
- Metatarsalgia •
- Pes planus •
- As a shock absorber •
- Keeping Achilles tendon intact •
- Wherever required •

### Size

XS

S

M

L

XL

XXL

24-25

26-27

28-29

30-31

32-33

34-35



## CORK INSOLE

OS 5509 Ç

It is manufactured, in compliance with the foot ergonomics, from epine cushioned cork and artificial leather with longitudinal arch and metatarsal support.

### Indication

- Fallen arch
- Pes planus
- Pes plano valgus



### Size

19

20

21

22

23

24

25

26

27

28

30

31

32

33

