



NELSONTYPE CERVICAL COLLAR **OS 7602**

Design with half-chin rest.
 Hypo-allergic 12mm. Made from plastazote material.
 Semi-rigidly shaped according to neck anatomy.
 It provides ventilaton for the skin with four holes on it.

Indication

- It is used
- In cervical disc deformations,
- In traumas
- In soft tissue injuries



Size

XS	S	M	L
25-28	29-33	34-38	39-46

OS 7603

NELSONTYPE CERVICAL COLLAR WITH CHIN REST

Design with full-chin rest. Hypo-allergic 12mm. Manufactured from 12 mm plastazote material.
 Semi-rigidly shaped according to the neck anatomy. It provides ventilaton for the skin with four holes on it.



Indication

- It is used
- In cervical disc deformations
- In traumas
- In soft tissue injuries
- In Osteo arthritis

Size

XS	S	M	L
25-28	29-33	34-38	39-46



SPONGE CERVICAL COLLAR **OS 7605**

It is manufactured from sponge material according to the cervical region anatomy.

Indication

- In slight neck injuries.



Size

S	M	L	XL
35-41	42-49	50-57	58-65



PHILADELPHIA TYPE CERVICAL COLLAR

OS 7601

It keeps the cervical region in the neutral position, preventing the hyperextension and flexion. It is shaped suitable to neck anatomy.

Indication

- In neck fractures
- In cervical disc deformations
- In trauma
- Emergency patient transfer after surgery



Size

XS	S	M	L
25-28	29-33	34-38	39-46

OS 7604

PHILADELPHIA TYPE CERVICAL COLLAR EXPOSING TRACHEA

It keeps the cervical region in the neutral position, preventing the hyperextension and flexion. It is shaped suitable to neck anatomy.



Indication

- In neck fractures
- In cervical disc deformations
- In trauma
- Emergency patient transfer after surgery

Size

S	M	L
29-33	34-38	39-46



SPONGE CERVICAL COLLAR WITH PLASTIC SUPPORT

OS 7605 A

It is manufactured from sponge material according to the cervical region anatomy.

Indication

- In slight neck injuries.



Size

S	M	L	XL
35-41	42-49	50-57	58-65



NELSON TYPE CERVICAL COLLAR WITH FABRIC COATING

OS 7602 A

Design with half-chin rest. Hypo-allergic 12mm. Made from plastazote material. Semi-rigidly shaped according to neck anatomy. It provides ventilation to the skin with four holes on it.

Indication

- In cervical disc deformations,
- In traumas
- In soft tissue injuries
- In Osteoarthritis.

Size

XS	S	M	L
25-28	29-33	34-38	39-46



OS 7608

HEIGHT ADJUSTABLE FIRST AID CERVICAL COLLAR

It is manufactured from rigid plastic according to the cervical region anatomy. Its height is adjustable three times. Its tracheal region is open for emergency intervention.



Indication

- Stabilization of the cervical region
- During emergency intervention

Size

XS	S	M	L
57-68	69-75	75-81	81-86



WIRE FRAME CERVICAL COLLAR

OS 7610

It is MR compatible. Extension level can be adjusted. The open front provides comfort and convenience for emergency intervention on.

Indication

- Traumatic cases
- After cervical surgery
- Cervical tendonitis
- Muscle spasm
- Disc hernia
- Osteoarthritis

Size

SMALL	LARGE		
25	33	34	46





CERVICAL COLLAR WITH CHIN REST

OS 7607

It is manufactured from rigid plastic according to the cervical region anatomy and height is adjustable. It ensures full stabilization with the chin rest

Indication

- In the neck hernias
- In trauma
- In post-operative stabilization cases



Size

S	M	L	XL	XXL
30-35	35-40	38-43	43-48	47-52

OS 7606

CERVICAL COLLAR WITHOUT CHIN REST

It is manufactured from rigid plastic according to the cervical region anatomy and height is adjustable.



Indication

- In the neck hernias
- In trauma
- In post-operative stabilization cases

Size

S	M	L	XL	XXL
30-35	35-40	38-43	43-48	47-52



CERVICAL COLLAR SLEEVE

OS 9170

It is made of cotton fabric to protect your skin and provide comfort. It is used over the necks. Manufactured from washable material.

Indication

- Used for comfort purposes



Size

STANDARD