



LUXURY ARMSLING

OS 2201

It is manufactured from laminated polyester sponge fabric the inner surface of which is cotton towel. Velcro fabric is used at the outer surface.. there is a belt to fix the sling to the waist.

Indication

- In case of forearm problems
- For resting purpose
- When the arm is required to be suspended
- In case of soft tissue injuries
- In cases where the arm must be immobilized.



Size

| | | | | |
|-------|-------|-------|-------|-------|
| XS | S | M | L | XL |
| 25-30 | 30-35 | 35-40 | 40-45 | 45-50 |

OS 2206

MESH ARMSLING

It is made of special perforated fabric. There is a sling which ensures that the arm is immobilized at the desired position. Also a child's size is available.



Indication

- When the arm is required to be suspended
- In soft tissue injuries
- When the arm needs to remain stationary

Size

| | | | | |
|-------|-------|-------|-------|-------|
| XS | S | M | L | XL |
| 25-30 | 30-35 | 35-40 | 40-45 | 45-50 |



CLOTH ARMSLING

OS 2209

It is made of black color fabric. Owing to the buckle on it, it is easy to use

Indication

- For resting purposes, in forearm problems
- When the arm is required to be suspended
- In soft tissue injuries
- In cases when the arm needs to remain motionless



Size

| | | | |
|-------|-------|-------|-------|
| S | M | L | XL |
| 30-32 | 34-36 | 38-40 | 40-42 |





ARMABDUCTION PILLOW 30°

OS 9002-A

It is manufactured from polyester sponge fabric the inner surface of which is cotton towel. It provides the immobilization at the 30° abduction of the shoulder. It Prevents post-operative adhesions of the tissues and helps mobilization and rehabilitation.

Indication

- Post-op impingement syndrome
- Rotator cuff tear repair
- Manipulation on the frozen shoulder under anesthesia
- Shoulder endoprosthesis



Size

STANDARD

OS 9002

ARMAADDUCTION PILLOW 45°-60°

It is manufactured from polyester sponge fabric the inner surface of which is cotton towel. It provides the immobilization at the 45° - 60° abduction of the shoulder. It Prevents post-operative adhesions of the tissues and helps mobilization and rehabilitation.



Indication

- Post-op impingement syndrome
- Rotator cuff tear repair
- Manipulation on the frozen shoulder under anesthesia
- Shoulder endoprosthesis

Size

STANDARD



SINGLE SHOULDER PAD

OS 1207

It is manufactured from neoprene fabric the inner surface of which is cotton. It supports shoulder.

Indication

- Shoulder subluxation
- In the dissipation of pain and edema
- For thermotherapy purposes



Size

S

M

L

XL

XXL

26-28

28-30

30-32

32-34

34-36



VELPEAUBANDAGE

OS 2208

It is manufactured from polyester sponge fabric the inner surface of which is cotton towel. Extra bandage was used to stabilize the shoulder joint,

Indication

- Shoulder subluxation fracture and subluxation
- For resting purposes in case of forearm problems.



Size

CHILD

ADULT

OS 9174

STIMSONBANDAGE

It is manufactured from a rigid fabric the inner of which is inelastic cotton fabric and the outer of which is velcro fabric. It is easy to use with the help of buckles and velcro on it.



Indication

- In the stabilization between Acampo-clavicle
- After Humerus operations



Size

STANDARD



CLAVICLE BANDAGE

OS 2205

It is manufactured from laminated polyester sponge fabric the inner surface of which is cotton. It has two plastic buckles. Its anatomical structure lets the clavicle to rest.

Indication

- Positioning the epiphyses in the collarbone fracture
- For stabilization purposes



Size

XS

S

M

L

XL

35-40

45-50

55-60

65-70

75-80

Size / Size

STANDAR

T UNISIZE