



## ELASTIC KNEE GUARD (NEOPRENE)

**OS 1402**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer

### Indication

- In painful condition of the knee
- Arthritis
- For Thermo therapy purposes



### Size

| S     | M     | L     | XL    | XXL   |
|-------|-------|-------|-------|-------|
| 21-23 | 23-25 | 25-27 | 27-29 | 29-31 |

**OS 1403**

## KNEE BRACE WITH PATELLASTABILIZER

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is a circular pad around patella to support it.



### Indication

- Arthritis
- Effusion
- In painful condition of the knee
- For Thermo therapy purposes

### Size

| S     | M     | L     | XL    | XXL   | XXXL  |
|-------|-------|-------|-------|-------|-------|
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



## KNEE BRACE WITH PATELLA AND LIGAMENT SUPPORT

**OS 1406**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is a circular pad around patella to support it. Mediolateral segments are supported by steel spiral.

### Indication

- Arthritis
- In cases of collateral ligament injuries
- In cross ligament injuries
- For Thermo therapy purposes



### Size

| S     | M     | L     | XL    | XXL   | XXXL  |
|-------|-------|-------|-------|-------|-------|
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



## KNEE BRACE WITH PATELLA AND LIGAMENT SUPPORT(WITH VELCRO)

**OS 1407**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is a circular pad around patella in order to support it. Mediolateral segments are supported by steel spiral. There are extra velcro at the upper calf and lower calf.

### Indication

- In Arthritis
- In cases of collateral ligament injuries
- In meniscus degenerations
- For Thermo therapy purposes

### Size

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| S     | M     | L     | XL    | XXL   | XXXL  |
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



**OS 1408**

## KNEE SUPPORT CROSSED LIGAMENTS

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is a circular pad around patella in order to support it. Mediolateral segments are supported by steel spiral. Extra strap was added with the aim of supporting cruciate ligaments



### Indication

- In anterior cruciate ligament injuries
- In collateral ligament injuries
- For thermotherapy purposes

### Size

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| S     | M     | L     | XL    | XXL   | XXXL  |
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



## KNEE SUPPORTCROSSED LIGAMENTS

**OS 1408 A**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is a circular pad around patella in order to support it. Mediolateral segments are supported by steel spiral. Extra strap was added with the aim of supporting cruciate ligaments. The knee brace is worn at the desired position by adjusting the velcro passing through the laminates on the knee brace.

### Indication

- In anterior cruciate ligament injuries
- In cases of collateral ligament injuries
- For Thermo therapy purposes

### Size

|       |       |       |       |       |
|-------|-------|-------|-------|-------|
| S     | M     | L     | XL    | XXL   |
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 |





## PATELLAR TENDONSTRAP

**OS 1404**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. It has an additional pad to stabilize the patellar tendonitis

### Indication

- Osgood schlatter
- Jumper's knee
- Tendonitis



### Size

STANDARD

**OS 1405**

## PATELLAR TENDONITIS KNEE BRACE

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. It has an additional pad to stabilize the tendonitis. The design that covers the knee completely, allows doing physical exercises.



### Indication

- Osgood schlatter
- Jumper's knee
- Tendonitis

### Size

| S     | M     | L     | XL    | XXL   |
|-------|-------|-------|-------|-------|
| 29-31 | 31-33 | 33-35 | 35-37 | 37-39 |



## STEEL HINGED KNEE BRACE WITH PATELLA STABILIZER

**OS 1409**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. It has an additional circular pad around patella to stabilize the patellar tendonitis. The mediolateral segments are supported by a steel hinge limited to 180°.

### Indication

- ACL-TCL
- In meniscus degenerations
- In collateral ligament injuries
- For thermotherapy purposes



### Size

| S     | M     | L     | XL    | XXL   | XXXL  |
|-------|-------|-------|-------|-------|-------|
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



## STOP HINGED KNEE BRACE WITH OPENABLE FRONT

**OS 1410**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. the front side can be opened completely. The mediolateral segments are supported by a steel hinge limited to 180°.



### Indication

- In anterior cruciate ligament injuries
- In cases of collateral ligament injuries

### Size

| S     | M     | L     | XL    | XXL   |
|-------|-------|-------|-------|-------|
| 33-36 | 36-39 | 39-42 | 42-45 | 45-48 |

**OS 1411**

## STEEL HINGED KNEE BRACE WITH CRUCIATE LIGAMENT SUPPORT

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. It has an additional circular pad around patella to stabilize the patellar tendonitis. The mediolateral segments are supported by a steel hinge limited to 180°. An additional strap is added in order to support cruciate ligaments.



### Indication

- In anterior cruciate ligament injuries
- In cases of collateral ligament injuries
- For Thermo therapy purposes

### Size

| S     | M     | L     | XL    | XXL   |
|-------|-------|-------|-------|-------|
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 |



## STOP HINGED KNEE BRACE WITH CROSSED STRAP AND LIGAMENT SUPPORT

**OS 1411 A**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. There is an additional circular pad around patella to stabilize the patellar tendonitis. he mediolateral segments are supported by a steel hinge limited to 180°. The knee brace is worn at the desired position by adjusting the velcro passing through the buckles on the knee brace.



### Indication

- In anterior cruciate ligament injuries
- In collateral ligament injuries
- For Thermo therapy purposes

### Size

| S     | M     | L     | XL    | XXL   | XXXL  |
|-------|-------|-------|-------|-------|-------|
| 30-32 | 32-34 | 35-37 | 37-39 | 40-42 | 43-45 |



## ADJUSTABLE ANGLE KNEE BRACE WITH PATELLA STABILIZER

**OS 1412**

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. The front side can completely open up. The mediolateral segments are supported by adjustable angle hinge.

### Indication

- In anterior cruciate ligament injuries
- In collateral ligament injuries
- During Post-op stage
- For Thermo therapy purposes



### Size

|          |          |          |           |            |             |
|----------|----------|----------|-----------|------------|-------------|
| <b>S</b> | <b>M</b> | <b>L</b> | <b>XL</b> | <b>XXL</b> | <b>XXXL</b> |
| 30-32    | 32-34    | 35-37    | 37-39     | 40-42      | 43-45       |

**OS 1414**

## ADJUSTABLE ANGLE KNEE BRACE WITH OPENABLE FRONT

It is manufactured from laminated neoprene fabric the inner surface of which is cotton towel layer. The mediolateral segments are supported by adjustable angle hinge. The front side can open fully.



### Indication

- In anterior cruciate ligament injuries
- In collateral ligament injuries
- During Post-op period
- For Thermo therapy purposes

### Size

|          |          |          |           |            |
|----------|----------|----------|-----------|------------|
| <b>S</b> | <b>M</b> | <b>L</b> | <b>XL</b> | <b>XXL</b> |
| 33-36    | 36-39    | 39-42    | 42-45     | 45-48      |



## KNEE BRACE WITH PADSUPPORT

**OS 1415**

It is manufactured from laminated neoprene fabric the inner surface of which is cotton towel layer. A pad is used at the patella region in order to absorb impacts. thanks to the velcro available at the upper and lower parts, the use of product is made easy.

### Indication

- To protect the knee cap from impacts



### Size

|          |          |          |           |            |
|----------|----------|----------|-----------|------------|
| <b>S</b> | <b>M</b> | <b>L</b> | <b>XL</b> | <b>XXL</b> |
| 30-32    | 32-34    | 35-37    | 37-39     | 40-42      |



## SHIN BRACE

OS 1413

It is manufactured from laminated neoprene fabric the inner surface of which is cotton towel layer. Velcro was added to facilitate the adaptation to lower calf.

### Indication

- Lower calf pain
- Lower calf tears
- Shin splint's
- For Thermo therapy purposes



### Size

|       |       |       |       |       |
|-------|-------|-------|-------|-------|
| S     | M     | L     | XL    | XXL   |
| 29-31 | 31-33 | 33-35 | 35-37 | 38-40 |

OS 1401

## THIGHBRACE

It is manufactured from laminated neoprene fabric the inner surface of which is cotton towel layer. It is designed to fit the thigh anatomy.



### Indication

- In the painful condition of the upper thigh
- In muscle rupture
- For Thermo therapy purposes



### Size

|       |       |       |       |       |
|-------|-------|-------|-------|-------|
| S     | M     | L     | XL    | XXL   |
| 39-41 | 42-44 | 45-47 | 48-50 | 51-54 |



## THERMAL TIGHT

OS 1301

It is manufactured, in compliance with the knee anatomy, from laminated neoprene fabric the inner surface of which is cotton towel layer. It relaxes the muscles by warming the surrounding area. The sweating feature allows cellulite to decrease and bases to shrink.

### Indication

- In cellulite treatment
- When exercising
- For Thermo therapy purposes



### Size

|       |       |       |       |        |         |
|-------|-------|-------|-------|--------|---------|
| S     | M     | L     | XL    | XXL    | XXXL    |
| 74-80 | 80-84 | 84-90 | 90-98 | 96-102 | 102-108 |