



NEOPRENE ANKLE GUARD (WITHOUT SUPPORT)

OS 1501

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel.

Indication

- In the Dissipation of post-distortion edema
- For thermotherapy purposes



Size

S	M	L	XL	XXL
21-22	22-23	23-24	24-25	25-26

OS 1503

ANKLE GUARD WITH LIGAMENT SUPPORT

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel. There are 4 pieces of steel spirals added in the mediolateral region in order to provide support for ligaments



Indication

- In cases where compression pad is required
- Application is required
- where stabilization is required
- In distortion

Size

S	M	L	XL	XXL
19-20	20-21	21-22	22-23	23-24



ANKLE GUARD WITH MALLEOLUS SUPPORT

OS 1505

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel. There are 4 pieces of steel spirals added in the mediolateral region in order to provide support for ligaments.

Indication

- In cases where compression pad application is required
- where stabilization is required
- In Distortion
- For thermotherapy purposes

Size

S	M	L	XL	XXL
19-20	20-21	21-22	22-23	23-24





NEOPRENE ANKLE GUARD (WITH ADDITIONAL BANDAGE)

OS 1502

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel. It has an extra ankle bandage.

Indication

- In cases where stabilization is required
- In the dissipation of post-distortion edema
- For thermo therapy purposes



Size

S	M	L	XL	XXL
19-20	20-21	21-22	22-23	23-24

OS 1504

ANKLE GUARD WITH PLASTIC SUPPORT

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel. Two plastic splints exposing the malleoli are added in order to provide support for ligaments in the mediolateral region.



Indication

- In mediolateral ligament injuries
- In cases where Stabilization is required
- In case of distortion
- For thermo therapy purposes

Size

S	M	L	XL	XXL
19-20	20-21	21-22	22-23	23-24



ANKLE GUARD WITH ACHILLES TENDONITIS SUPPORT

OS 1506

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel. 4 pieces of steel spirals are added for supporting the ligaments.

Indication

- In cases where compression pad application is required
- In cases where stabilization is required
- In case of distortion
- For thermo therapy purposes



Size

S	M	L	XL	XXL
19-20	20-21	21-22	22-23	23-24





ANKLE STABILIZER (WITH SPONGE PAD)

OS 1510

It is manufactured, in compliance with the ankle anatomy, from laminated neoprene fabric the inner surface of which is cotton towel.

Indication

- When the knee is painful
- Arthritis
- For thermo therapy purposes



Size

STANDARD

OS 1508

DORSIFLEXION STRAP OS1508-A SHOE APPARATUS

It is manufactured anatomically from two separate pieces such as hardened velcro fabric and buckle in order to prevent Drop of the ankle. Use apparatus for shoe is available as an option.



Indication

In case of foot drop syndromes •

Size

XS-S

M-L

XL-XXL

25-30

30-35

35-40



ANKLE BANDAGE

OS 1509

It is manufactured, in compliance with ankle anatomy from cotton fabric.

Indication

- In the dissipation of post-distortion edema
- For thermo therapy purposes
- In the slight stabilization of ankle



Size

STANDARD



ANKLE STABILIZER (AIR PAD)

OS 1515

It is manufactured from polyethylene material inside of which there is an inflatable air pad.

Indication

- Post-op ligament injuries
- Post-sprains
- In sports injuries
- Post-op plaster cast stabilization
- Inversion and eversion injuries

Size

STANDARD



OS 1513

ANKLE STABILIZATION ORTHESIS

It is manufactured from polyethylene material the inside of which is cotton fabric and the outer surface of which is supported with velcro. It offers easy and practical use owing to its anatomical structure.



Indication

- Post-op ligament injuries
- post-sprains
- In sports injuries
- Post-op plaster cast stabilization
- Inversion and eversion injuries

Size

STANDARD



HALLUX VALGUS NIGHT SPLINT

OS 1507

It is manufactured as right or left to separate the big toe from the other toes.

Indication

- Hallux Valgus disorder

Size

STANDARD



< LEFT RIGHT >